

Surface Disinfectants

Perimo RTU

All Surface Disinfectant/Cleaner

Perimo RTU™ uses quaternary ammonium to effectively kill Tuberculosis, HIV, Hepatitis B and C, Polio and drug resistant staph (VRE/MRSA). Perimo RTU does not absorb through the skin and will not stain, damage instruments or hard surfaces. It can be used for immersion, as a surface spray or in an ultrasonic unit.

- 5 minute TB kill time, 3 minutes for most major viruses including HBV, HCV, MRSA, Rhinovirus, Rotavirus, 1 minute for HIV
 - Ready to use Perimo is easy to use and does not require mixing or measuring
 - Non-staining and non-flammable
 - 2 year shelf life
 - Versatile - can be used to soak instruments or in an ultrasonic cleaner
 - Eliminates toxicity, inhalation and skin absorption dangers of phenol
 - Comes in ready to use quarts, gallons and five gallon container.
-
- B025 - Quart
 - B026 - Gallon
 - B027 - 5 Gallon



**For Ordering Information
Contact CHARTECH INC.
1-877-658-1849
Fax 770-442-2843
sales@chartech.net**